



Sterile Container System & Baskets



We produce the best for you

Classic & Bio-Barrier Model



OLUŞUM[®]

OLUŞUM CERRAHİ ALETLER SAN. VE TİC. LTD. ŞTİ.
OLUŞUM SURGICAL INSTRUMENTS INDUSTRY AND TRADE CO.
Aşağıçinik Beldesi Aşmek Sanayi Sitesi No: 5-6 Tekkeköy/SAMSUN
Tel : 90 (362) 264 81 97 - 264 83 17 Fax : 90 (362) 264 84 27
www.olusumcerrahi.com E-mail : bilgi@olusumcerrahi.com

ISO 13485:2003

CE




ISO 9001 : 2008




- According to the production programme Oluşum's rights is reserved in production to improve, modify or delete the the products declared in this catalog.

- The weight information given in our catalog is approximate value.



The holding pins and plates for separation wire baskets or desire



The silicon mat, for placing the micro instruments.



Silicone handle has available in six different colors and detachable by user.




The drape retainer.



Spiral holder




Wire bottom, balanced placement and steam circulation safety.




The fixation clamps for the cylindrical instruments in different sizes.




The method of mounting - de-mounting the safety lid of the fix model container.




1. Before sterilization the valves are in relax position.




2. During the sterilization by effect of vacuum the valves are opened completely and the dry air goes out of the container.



3. At the end the sterilization cycle the valves are opened in the reverse way by effect of vacuum and the sterilization cycle is completed.



4. At the end of the sterilization cycle the sterilizer equates the pressure level and the valves of the Barrier Model Container are closed.



5. After sterilization the valves are in relax position.









Sterile Container System




1/1 Full Size

Classic Model


Specifications Sizes	Lid (Perforated) 	Lid (Perforated) 	Safety Lid 	Safety Lid 
	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)	Box (Perforated)
580 X 280 X 100 mm 	30-010  30-011  30-012  30-013  30-014  30-015 	30-040  30-041  30-042  30-043  30-044  30-045 	30-150  30-151  30-152  30-153  30-154  30-155 	30-180  30-181  30-182  30-183  30-184  30-185 
580 X 280 X 135 mm 	30-016  30-017  30-018  30-019  30-020  30-021 	30-046  30-047  30-048  30-049  30-050  30-051 	30-156  30-157  30-158  30-159  30-160  30-161 	30-186  30-187  30-188  30-189  30-190  30-191 
580 X 280 X 150 mm 	30-022  30-023  30-024  30-025  30-026  30-027 	30-052  30-053  30-054  30-055  30-056  30-057 	30-162  30-163  30-164  30-165  30-166  30-167 	30-192  30-193  30-194  30-195  30-196  30-197 
580 X 280 X 200 mm 	30-028  30-029  30-030  30-031  30-032  30-033 	30-058  30-059  30-060  30-061  30-062  30-063 	30-168  30-169  30-170  30-171  30-172  30-173 	30-198  30-199  30-200  30-201  30-202  30-203 
580 X 280 X 260 mm 	30-034  30-035  30-036  30-037  30-038  30-039 	30-064  30-065  30-066  30-067  30-068  30-069 	30-174  30-175  30-176  30-177  30-178  30-179 	30-204  30-205  30-206  30-207  30-208  30-209 
Ps: All container measurements cover the outer size of the lids.				
				



CODE	OUTSIDE	INSIDE	WEIGHT
30-100	595 X 275 X 95 mm	543 X 260 X 75 mm	2.15 kg
30-101	595 X 275 X 116 mm	543 X 260 X 96 mm	2.30 kg
30-102	595 X 275 X 135 mm	543 X 260 X 115 mm	2.40 kg
30-103	595 X 275 X 185 mm	543 X 260 X 165 mm	2.75 kg
30-104	595 X 275 X 243 mm	543 X 260 X 223 mm	3.25 kg



CODE	OUTSIDE	INSIDE	WEIGHT
30-105	595 X 275 X 95 mm	543 X 260 X 90 mm	1.70 kg
30-106	595 X 275 X 116 mm	543 X 260 X 111 mm	2.05 kg
30-107	595 X 275 X 135 mm	543 X 260 X 130 mm	2.15 kg
30-108	595 X 275 X 185 mm	543 X 260 X 180 mm	2.60 kg
30-109	595 X 275 X 243 mm	543 X 260 X 238 mm	3.10 kg

**Perforated Lids****Safety Lids****Removal Lids (Non Perforated)**

OLUŞUM[®]







Sterile Container System




3/4 Quarter Size

Classic Model

Specifications Sizes	Lid (Perforated)	Lid (Perforated)	Safety Lid	Safety Lid
	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)	Box (Perforated)
465 X 280 X 100 mm 	30-240  30-241  30-242  30-243  30-244  30-245 	30-260  30-261  30-262  30-263  30-264  30-265 	30-350  30-351  30-352  30-353  30-354  30-355 	30-370  30-371  30-372  30-373  30-374  30-375 
465 X 280 X 135 mm 	30-246  30-247  30-248  30-249  30-250  30-251 	30-266  30-267  30-268  30-269  30-270  30-271 	30-356  30-357  30-358  30-359  30-360  30-361 	30-376  30-377  30-378  30-379  30-380  30-381 
465 X 280 X 150 mm 	30-252  30-253  30-254  30-255  30-256  30-257 	30-272  30-273  30-274  30-275  30-276  30-277 	30-362  30-363  30-364  30-365  30-366  30-367 	30-382  30-383  30-384  30-385  30-386  30-387 
Ps: All container measurements cover the outer size of the lids.				
				



CODE	OUTSIDE	INSIDE	WEIGHT
30-300	475 X 275 X 95 mm	422 X 260 X 90 mm	1.65 kg
30-301	475 X 275 X 116 mm	422 X 260 X 96 mm	1.80 kg
30-302	475 X 275 X 137 mm	422 X 260 X 117 mm	1.95 kg




CODE	OUTSIDE	INSIDE	WEIGHT
30-305	475 X 275 X 95 mm	422 X 260 X 90 mm	1.50 kg
30-306	475 X 275 X 116 mm	422 X 260 X 111 mm	1.65 kg
30-307	475 X 275 X 137 mm	422 X 260 X 132 mm	1.75 kg




Perforated Lids



Safety Lids




Removal Lids (Non Perforated)



GREY
YELLOW
GREEN
BLUE
RED
BLACK









Sterile Container System




1/2 Half Size

Classic Model

Specifications Sizes	Lid (Perforated)	Lid (Perforated)	Safety Lid	Safety Lid
	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)	Box (Perforated)
285 X 280 X 100 mm 	30-500  30-501  30-502  30-503  30-504  30-505 	30-530  30-531  30-532  30-533  30-534  30-535 	30-640  30-641  30-642  30-643  30-644  30-645 	30-670  30-671  30-672  30-673  30-674  30-675 
285 X 280 X 135 mm 	30-506  30-507  30-508  30-509  30-510  30-511 	30-536  30-537  30-538  30-539  30-540  30-541 	30-646  30-647  30-648  30-649  30-650  30-651 	30-676  30-677  30-678  30-679  30-680  30-681 
285 X 280 X 150 mm 	30-512  30-513  30-514  30-515  30-516  30-517 	30-542  30-543  30-544  30-545  30-546  30-547 	30-652  30-653  30-654  30-655  30-656  30-657 	30-682  30-683  30-684  30-685  30-686  30-687 
285 X 280 X 200 mm 	30-518  30-519  30-520  30-521  30-522  30-523 	30-548  30-549  30-550  30-551  30-552  30-553 	30-658  30-659  30-660  30-661  30-662  30-663 	30-688  30-689  30-690  30-691  30-692  30-693 
285 X 280 X 260 mm 	30-524  30-525  30-526  30-527  30-528  30-529 	30-554  30-555  30-556  30-557  30-558  30-559 	30-664  30-665  30-666  30-667  30-668  30-669 	30-694  30-695  30-696  30-697  30-698  30-699 
Ps: All container measurements cover the outer size of the lids.				
   				
   				



CODE	OUTSIDE	INSIDE	WEIGHT
30-600	306 X 272 X 95 mm	253 X 258 X 75 mm	1.25 kg
30-601	306 X 272 X 114 mm	253 X 258 X 94 mm	1.40 kg
30-602	306 X 272 X 135 mm	253 X 258 X 115 mm	1.55 kg
30-603	306 X 272 X 184 mm	253 X 258 X 164 mm	1.80 kg
30-604	306 X 272 X 241 mm	253 X 258 X 221 mm	2.10 kg



CODE	OUTSIDE	INSIDE	WEIGHT
30-605	306 X 272 X 95 mm	253 X 258 X 90 mm	1.15 kg
30-606	306 X 272 X 114 mm	253 X 258 X 109 mm	1.25 kg
30-607	306 X 272 X 135 mm	253 X 258 X 130 mm	1.35 kg
30-608	306 X 272 X 184 mm	253 X 258 X 179 mm	1.55 kg
30-609	306 X 272 X 241 mm	253 X 258 X 236 mm	2.00 kg



Perforated Lids



30-630

30-631


30-632

30-633

30-634

30-635

Safety Lids



30-700

30-701


30-702

30-703

30-704

30-705

Removal Lids (Non Perforated)



30-620

30-621

30-622

30-623

30-624

30-625



GREY

YELLOW

GREEN

BLUE

RED

BLACK

OLUŞUM[®]






Sterile Container System




XL Size

Classic Model

Specifications Sizes	Lid (Perforated)	Lid (Perforated)
	Box (Non-Perforated)	Box (Perforated)
590 X 370 X 120 mm 	30-136  30-137  30-138  30-139  30-140  30-141 	30-210  30-211  30-212  30-213  30-214  30-215 
590 X 370 X 180 mm 	30-142  30-143  30-144  30-145  30-146  30-147 	30-216  30-217  30-218  30-219  30-220  30-221 
Ps: All container measurements cover the outer size of the lids.		
		



CODE	OUTSIDE	INSIDE	WEIGHT
30-400	595 X 365 X 105 mm	530 X 350 X 85 mm	2.50 kg
30-401	595 X 365 X 167 mm	530 X 350 X 145 mm	3.25 kg



CODE	OUTSIDE	INSIDE	WEIGHT
30-402	595 X 365 X 105 mm	543 X 350 X 100 mm	2.20 kg
30-403	595 X 365 X 167 mm	530 X 350 X 162 mm	3.00 kg



Perforated Lids



- 30-326**
- 30-327**
- 30-328**
- 30-329**
- 30-330**
- 30-331**




- GREY
- YELLOW
- GREEN
- BLUE
- RED
- BLACK

OLUŞUM[®]








Sterile Container System




1/1 Full Size

Bio-Barrier Model

Specifications		Lid (with barrier)	Lid (with barrier)	Safety Lid Lid (with barrier)	Safety Lid Lid (with barrier)
Sizes		 Box (Non-Perforated)	 Box (with barrier)	 Box (Non-Perforated)	 Box (with barrier)
580 X 280 X h mm 	31-010  31-011  31-012  31-013  31-014  31-015  h= 115 mm	31-040  31-041  31-042  31-043  31-044  31-045  h= 135 mm	31-150  31-151  31-152  31-153  31-154  31-155  h= 128 mm	31-180  31-181  31-182  31-183  31-184  31-185  h= 145 mm	
580 X 280 X h mm 	31-016  31-017  31-018  31-019  31-020  31-021  h= 140 mm	31-046  31-047  31-048  31-049  31-050  31-051  h= 163 mm	31-156  31-157  31-158  31-159  31-160  31-161  h= 153 mm	31-186  31-187  31-188  31-189  31-190  31-191  h= 170 mm	
580 X 280 X h mm 	31-022  31-023  31-024  31-025  31-026  31-027  h= 155 mm	31-052  31-053  31-054  31-055  31-056  31-057  h= 175 mm	31-162  31-163  31-164  31-165  31-166  31-167  h= 168 mm	31-192  31-193  31-194  31-195  31-196  31-197  h= 185 mm	
580 X 280 X h mm 	31-028  31-029  31-030  31-031  31-032  31-033  h= 205 mm	31-058  31-059  31-060  31-061  31-062  31-063  h= 235 mm	31-168  31-169  31-170  31-171  31-172  31-173  h= 218 mm	31-198  31-199  31-200  31-201  31-202  31-203  h= 245 mm	
580 X 280 X h mm 	31-034  31-035  31-036  31-037  31-038  31-039  h= 265 mm	31-064  31-065  31-066  31-067  31-068  31-069  h= 295 mm	31-174  31-175  31-176  31-177  31-178  31-179  h= 278 mm	31-204  31-205  31-206  31-207  31-208  31-209  h= 305 mm	
Ps: All container measurements cover the outer size of the lids.					
					



CODE	OUTSIDE	INSIDE	WEIGHT
31-100	595 X 275 X 115 mm	543 X 260 X 90 mm	2.70 kg
31-101	595 X 275 X 136 mm	543 X 260 X 111 mm	3.20 kg
31-102	595 X 275 X 155 mm	543 X 260 X 130 mm	3.40 kg
31-103	595 X 275 X 205 mm	543 X 260 X 180 mm	4.10 kg
31-104	595 X 275 X 263 mm	543 X 260 X 238 mm	4.90 kg



CODE	OUTSIDE	INSIDE	WEIGHT
31-105	595 X 275 X 95 mm	543 X 260 X 90 mm	1.70 kg
31-106	595 X 275 X 116 mm	543 X 260 X 111 mm	2.05 kg
31-107	595 X 275 X 135 mm	543 X 260 X 130 mm	2.15 kg
31-108	595 X 275 X 185 mm	543 X 260 X 180 mm	2.60 kg
31-109	595 X 275 X 243 mm	543 X 260 X 238 mm	3.10 kg


**Lid Barrier system****Safety Lids**

GREY
YELLOW
GREEN
BLUE
RED
BLACK

OLUŞUM[®]








Sterile Container System




3/4 Quarter Size

Bio-Barrier Model

Specifications Sizes	Lid (with barrier)  Box (Non-Perforated)	Lid (with barrier)  Box (with barrier)	Safety Lid Lid (with barrier)  Box (Non-Perforated)	Safety Lid Lid (with barrier)  Box (with barrier)
465 X 280 X h mm 	31-240  31-241  31-242  31-243  31-244  31-245  h= 115 mm	31-260  31-261  31-262  31-263  31-264  31-265  h= 135 mm	31-350  31-351  31-352  31-353  31-354  31-355  h= 128 mm	31-370  31-371  31-372  31-373  31-374  31-375  h= 150 mm
465 X 280 X h mm 	31-246  31-247  31-248  31-249  31-250  31-251  h= 140 mm	31-266  31-267  31-268  31-269  31-270  31-271  h= 160 mm	31-356  31-357  31-358  31-359  31-360  31-361  h= 153 mm	31-376  31-377  31-378  31-379  31-380  31-381  h= 175 mm
465 X 280 X h mm 	31-252  31-253  31-254  31-255  31-256  31-257  h= 155 mm	31-272  31-273  31-274  31-275  31-276  31-277  h= 175 mm	31-362  31-363  31-364  31-365  31-366  31-367  h= 168 mm	31-382  31-383  31-384  31-385  31-386  31-387  h= 185 mm
Ps: All container measurements cover the outer size of the lids.				
   				
   				




CODE	OUTSIDE	INSIDE	WEIGHT
31-300	475 X 275 X 115 mm	422 X 260 X 90 mm	1.95 kg
31-301	475 X 275 X 136 mm	422 X 260 X 111 mm	2.10 kg
31-302	475 X 275 X 157 mm	422 X 260 X 132 mm	2.20 kg




CODE	OUTSIDE	INSIDE	WEIGHT
31-305	475 X 275 X 95 mm	422 X 260 X 90 mm	1.50 kg
31-306	475 X 275 X 116 mm	422 X 260 X 111 mm	1.65 kg
31-307	475 X 275 X 137 mm	422 X 260 X 132 mm	1.75 kg



Lid Barrier System




Safety Lids



OLUŞUM®







Sterile Container System




1/2 Half Size

Bio-Barrier Model

Specifications Sizes	Lid (with barrier)		Lid (with barrier)		Safety Lid Lid (with barrier)		Safety Lid Lid (with barrier)	
	 Box (Non-Perforated)	 Box (with barrier)	 Box (Non-Perforated)	 Box (with barrier)				
285 X 280 X h mm 	31-500  31-501  31-502  31-503  31-504  31-505 	h= 105 mm	31-530  31-531  31-532  31-533  31-534  31-535 	h= 135 mm	31-640  31-641  31-642  31-643  31-644  31-645 	h= 128 mm	31-670  31-671  31-672  31-673  31-674  31-675 	h= 145 mm
285 X 280 X h mm 	31-506  31-507  31-508  31-509  31-510  31-511 	h= 140 mm	31-536  31-537  31-538  31-539  31-540  31-541 	h= 160 mm	31-646  31-647  31-648  31-649  31-650  31-651 	h= 153 mm	31-676  31-677  31-678  31-679  31-680  31-681 	h= 170 mm
285 X 280 X h mm 	31-512  31-513  31-514  31-515  31-516  31-517 	h= 155 mm	31-542  31-543  31-544  31-545  31-546  31-547 	h= 175 mm	31-652  31-653  31-654  31-655  31-656  31-657 	h= 168 mm	31-682  31-683  31-684  31-685  31-686  31-687 	h= 185 mm
285 X 280 X h mm 	31-518  31-519  31-520  31-521  31-522  31-523 	h= 205 mm	31-548  31-549  31-550  31-551  31-552  31-553 	h= 235 mm	31-658  31-659  31-660  31-661  31-662  31-663 	h= 218 mm	31-688  31-689  31-690  31-691  31-692  31-693 	h= 245 mm
285 X 280 X h mm 	31-524  31-525  31-526  31-527  31-528  31-529 	h= 265 mm	31-554  31-555  31-556  31-557  31-558  31-559 	h= 295 mm	31-664  31-665  31-666  31-667  31-668  31-669 	h= 278 mm	31-694  31-695  31-696  31-697  31-698  31-699 	h= 305 mm
Ps: All container measurements cover the outer size of the lids.								
   								
   								




CODE	OUTSIDE	INSIDE	WEIGHT
31-600	306 X 272 X 115 mm	253 X 258 X 90 mm	1.60 kg
31-601	306 X 272 X 134 mm	253 X 258 X 109 mm	1.70 kg
31-602	306 X 272 X 155 mm	253 X 258 X 130 mm	1.80 kg
31-603	306 X 272 X 204 mm	253 X 258 X 179 mm	2.00 kg
31-604	306 X 272 X 261 mm	253 X 258 X 236 mm	2.45 kg



CODE	OUTSIDE	INSIDE	WEIGHT
31-605	306 X 272 X 95 mm	253 X 258 X 90 mm	1.15 kg
31-606	306 X 272 X 114 mm	253 X 258 X 109 mm	1.25 kg
31-607	306 X 272 X 135 mm	253 X 258 X 130 mm	1.35 kg
31-608	306 X 272 X 184 mm	253 X 258 X 179 mm	1.55 kg
31-609	306 X 272 X 241 mm	253 X 258 X 236 mm	2.00 kg




Lid Barrier system



31-630	
31-631	
31-632	
31-633	
31-634	
31-635	

Safety Lids



31-412	
31-413	
31-414	
31-415	
31-416	
31-417	




GREY
YELLOW
GREEN
BLUE
RED
BLACK

OLUŞUM[®]






Sterile Container System




XL Size

Bio-Barrier Model

Specifications Sizes	Lid (with barrier) Box (Non-Perforated)		Lid (with barrier) Box (with barrier)	
	590 X 370 X h mm 		590 X 370 X h mm 	
	31-711  31-712  31-713  31-714  31-715  31-716 	h = 125 mm	31-741  31-742  31-743  31-744  31-745  31-746 	h = 135 mm
	31-717  31-718  31-719  31-720  31-721  31-722 	h = 185 mm	31-747  31-748  31-749  31-750  31-751  31-752 	h = 195 mm
Ps: All container measurements cover the outer size of the lids.				
 				
 				



CODE	OUTSIDE	INSIDE	WEIGHT
31-700	595 X 365 X 125 mm	530 X 350 X 100 mm	2.620 kg
31-701	595 X 365 X 187 mm	530 X 350 X 162 mm	3.440 kg



CODE	OUTSIDE	INSIDE	WEIGHT
31-702	595 X 365 X 105 mm	530 X 350 X 100 mm	2.20 kg
31-703	595 X 365 X 167 mm	530 X 350 X 167 mm	3.00 kg




Lid Barrier system




- 31-771**
- 31-772**
- 31-773**
- 31-774**
- 31-775**
- 31-776**




- GREY
- YELLOW
- GREEN
- BLUE
- RED
- BLACK




1. Before sterilization the valves are in relax position.




2. During the sterilization by effect of vacuum the valves are opened completely and the dry air goes out of the container.



3. At the end the sterilization cycle the valves are opened in the reverse way by effect of vacuum and the sterilization cycle is completed.




4. At the end of the sterilization cycle the sterilizer equates the pressure level and the valves of the Barrier Model Container are closed.



5. After sterilization the valves are in relax position.



Sterile Container System




Dental Mini & Special Model

Specifications Sizes	Lid (Perforated)	Lid (Perforated)	Lid (Non Perforated)		
	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)
320 X 190 X 50 mm 	30-1800 30-1801 30-1802 30-1803 30-1804 30-1805	30-1820 30-1821 30-1822 30-1823 30-1824 30-1825	30-1840 30-1841 30-1842 30-1843 30-1844 30-1845	30-1860	30-1865
320 X 190 X 65 mm 	30-1806 30-1807 30-1808 30-1809 30-1810 30-1811	30-1826 30-1827 30-1828 30-1829 30-1830 30-1831	30-1846 30-1847 30-1848 30-1849 30-1850 30-1851	30-1861	30-1866
320 X 190 X 130 mm 	30-1812 30-1813 30-1814 30-1815 30-1816 30-1817	30-1832 30-1833 30-1834 30-1835 30-1836 30-1837	30-1852 30-1853 30-1854 30-1855 30-1856 30-1857	30-1862	30-1867

Ps: All container measurements cover
the outer size of the lids.




33-1070 274 x 172 x 43 mm
WIRE BASKET



33-611 274 x 172 x 43 mm
WIRE BASKET WITH LID.



33-609 273 x 180 x 13 mm
STORAGE FRAME
for dental model




32-529 260 x 180 mm
SILICONE MAT








33-610 277 x 175 mm
DRAPE RETAINER



Perforated Bottoms



Solid Bottoms



Perforated Lids



33-615


Single use paper filter 100 pcs.
With indicator No marking.




32-743

Teflon filter
Number of sterilization cycle > 1200

Bio-Barrier Lids



Removal Lids



Specifications Sizes	Lid (Perforated)	Lid (Perforated)	Lid (Non Perforated)		
	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)
315 X 135 X 40 mm 	30-1900 30-1901 30-1902 30-1903 30-1904 30-1905	30-1920 30-1921 30-1922 30-1923 30-1924 30-1925	30-1940 30-1941 30-1942 30-1943 30-1944 30-1945	30-1960	30-1965
315 X 135 X 70 mm 	30-1906 30-1907 30-1908 30-1909 30-1910 30-1911	30-1926 30-1927 30-1928 30-1929 30-1930 30-1931	30-1946 30-1947 30-1948 30-1949 30-1950 30-1951	30-1961	30-1966
315 X 135 X 100 mm 	30-1912 30-1913 30-1914 30-1915 30-1916 30-1917	30-1932 30-1933 30-1934 30-1935 30-1936 30-1937	30-1952 30-1953 30-1954 30-1955 30-1956 30-1957	30-1962	30-1967

Ps: All container measurements cover
the outer size of the lids.



33-1073 267 x 125 x 50 mm
WIRE BASKET



32-993 267 x 125 x 50 mm
WIRE BASKET WITH LID.



32-990 270 x 125 mm
DRAPE RETAINER



32-992 270 x 125 mm
WIRE BASE


32-991

Portable carrier for small containers.
Wire suspenders adjustable as required.
Adjusting key of suspender is included.
Width: 325 mm, Height: 328 mm, Depth: 145 mm








32-520 275 x 125 mm
SILICONE MAT

Perforated Bottoms



Solid Bottoms


**32-761**

Single use paper filter 100 pcs.
With indicator. No marking.


**32-744**

Teflon filter
Number of sterilization scycle > 1200

Perforated Lid



Removal Lids



IMPLANT Model Containers

Specifications / Sizes	Lid (Perforated)	Lid (Perforated)	Lid (Non Perforated)	NON FILTER
	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)	
300 X 149 X 65 mm 	29-010  29-011  29-012  29-013  29-014  29-015 	29-022  29-023  29-024  29-025  29-026  29-027 	29-034  29-035  29-036  29-037  29-038  29-039 	
500 X 169 X 75 mm 	29-016  29-017  29-018  29-019  29-020  29-021 	29-028  29-029  29-030  29-031  29-032  29-033 	29-040  29-041  29-042  29-043  29-044  29-045 	

Ps: All container measurements cover the outer size of the lids.

MINI IMPLANT Model Containers

Specifications / Sizes	Lid (Perforated)	Lid (Perforated)	Lid (Non Perforated)	NON FILTER
	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)	
160 X 70 X 60 mm 	29-200  29-201  29-202  29-203  29-204  29-205 	29-212  29-213  29-214  29-215  29-216  29-217 	29-224  29-225  29-226  29-227  29-228  29-229 	
200 X 145 X 60 mm 	29-206  29-207  29-208  29-209  29-210  29-211 	29-218  29-219  29-220  29-221  29-222  29-223 	29-230  29-231  29-232  29-233  29-234  29-235 	

Ps: All container measurements cover the outer size of the lids.

**BRUSH Model
Containers**





-  **32-044** Black
-  **32-045** Red
-  **32-046** Blue
-  **32-047** Green
-  **32-048** Yellow
-  **32-049** Grey

Container for brush 280 x 285 x 85 mm with wall mounting fixtures, Maximum 20 each brushes can be stored container. To take the brush, please move the handle to right or left direction under the container.

non filter

Specifications / Sizes	Lid (Perforated)	Lid (Perforated)
	Box (Non-Perforated)	Box (Perforated)
250 X 60 X 50 mm 	29-060 29-061 29-062 29-063 29-064 29-065	29-072 29-073 29-074 29-075 29-076 29-077
450 X 70 X 70 mm 	29-066 29-067 29-068 29-069 29-070 29-071	29-078 29-079 29-080 29-081 29-082 29-083

Ps: All container measurements cover the outer size of the lids.



SCOPY Model Containers

Specifications / Sizes	Lid (Perforated)	Lid (Perforated)
	Box (Non-Perforated)	Box (Perforated)
250 X 60 X 50 mm 	29-100 29-101 29-102 29-103 29-104 29-105	29-115 29-116 29-117 29-118 29-119 29-120
450 X 70 X 70 mm 	29-106 29-107 29-108 29-109 29-110 29-111	29-121 29-122 29-123 29-124 29-125 29-126

Ps: All container measurements cover the outer size of the lids.



SCOPY Model Containers
(with silicon stand, lock, label)

Specifications Sizes	Lid (Perforated)	Lid (Perforated)	Lid (Non Perforated)		
	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)
285 X 280 X 55 mm 	30-1710  30-1711  30-1712  30-1713  30-1714  30-1715 	30-1730  30-1731  30-1732  30-1733  30-1734  30-1735 	30-1750  30-1751  30-1752  30-1753  30-1754  30-1755 	30-1770 	30-1775 
285 X 280 X 85 mm 	30-1716  30-1717  30-1718  30-1719  30-1720  30-1721 	30-1736  30-1737  30-1738  30-1739  30-1740  30-1741 	30-1756  30-1757  30-1758  30-1759  30-1760  30-1761 	30-1771 	30-1776 

Ps: All container measurements cover
the outer size of the lids.




32-490 255 x 250 x 30 mm
32-491 255 x 250 x 50 mm
 WIRE BASKET




32-482 255 x 255 x 30 mm
 WIRE BASE



32-480 255 x 255 mm
 DRAPE RETAINER



32-481 255 x 255 x 30 mm
 DRAPE RETAINER



32-521 255 x 255 mm
 SILICON MAT




32-740
 Single use paper filter 100 pcs.




32-742
 Teflon Filter
 Number of sterilization cycle > 1200


Perforated Bottoms



Solid Bottoms




Perforated Lids



Bio-Barrier Lids




Safety Lids



Bio-Barrier Safety Lids




Removal Lids



OLUŞUM[®]



Sterile Container System



Baskets



33-1010	540 x 250 x 30 mm
33-1011	540 x 250 x 50 mm
33-1012	540 x 250 x 70 mm
33-1013	540 x 250 x 100 mm

wire basket with lid

33-1014 ONLY LID (for baskets)



32-460	540 x 250 x 30 mm
32-461	540 x 250 x 50 mm
32-462	540 x 250 x 70 mm
32-463	540 x 250 x 100 mm

wire basket



33-1020	480 x 250 x 30 mm
33-1021	480 x 250 x 50 mm
33-1022	480 x 250 x 70 mm
33-1023	480 x 250 x 100 mm

wire basket with lid

33-1024 ONLY LID (for baskets)



32-450	480 x 250 x 30 mm
32-451	480 x 250 x 50 mm
32-452	480 x 250 x 70 mm
32-453	480 x 250 x 100 mm

wire basket



33-1030	405 x 250 x 30 mm
33-1031	405 x 250 x 50 mm
33-1032	405 x 250 x 70 mm
33-1033	405 x 250 x 100 mm

wire basket with lid

33-1034 ONLY LID (for baskets)



32-1500	405 x 250 x 30 mm
32-1501	405 x 250 x 50 mm
32-1502	405 x 250 x 70 mm
32-1503	405 x 250 x 100 mm

wire basket



33-1040	255 x 250 x 30 mm
33-1041	255 x 250 x 50 mm
33-1042	255 x 250 x 70 mm
33-1043	255 x 250 x 100 mm

wire basket with lid

33-1045 ONLY LID (for baskets)



32-490	255 x 250 x 30 mm
32-491	255 x 250 x 50 mm
32-492	255 x 250 x 70 mm
32-493	255 x 250 x 100 mm

wire basket

Wire Basket With Perforated Metal Plate





WIRE BASKET

32-461 P	540 x 250 x 50 mm
32-462 P	540 x 250 x 70 mm
32-463 P	540 x 250 x 100 mm
32-451 P	480 x 250 x 50 mm
32-452 P	480 x 250 x 70 mm
32-453 P	480 x 250 x 100 mm



WIRE BASKET WITH LID

33-1011 P	540 x 250 x 50 mm
33-1012 P	540 x 250 x 70 mm
33-1013 P	540 x 250 x 100 mm
33-1014 P	ONLY LID (for baskets)
33-1021 P	480 x 250 x 50 mm
33-1022 P	480 x 250 x 70 mm
33-1023 P	480 x 250 x 100 mm
33-1024 P	ONLY LID (for baskets)




WIRE BASKET

32-1501 P	405 x 250 x 50 mm
32-1502 P	405 x 250 x 70 mm
32-1503 P	405 x 250 x 100 mm



WIRE BASKET WITH LID

33-1031 P	405 x 250 x 50 mm
33-1032 P	405 x 250 x 70 mm
33-1033 P	405 x 250 x 100 mm
33-1034 P	ONLY LID (for baskets)




WIRE BASKET


32-491 P	255 x 250 x 50 mm
32-492 P	255 x 250 x 70 mm
32-493 P	255 x 250 x 100 mm




WIRE BASKET WITH LID

33-1041 P	255 x 250 x 50 mm
33-1042 P	255 x 250 x 70 mm
33-1043 P	255 x 250 x 100 mm
33-1045 P	ONLY LID (for baskets)



**33-1050** 670 x 80 x 52 mm**33-1051** 460 x 80 x 52 mm**33-1052** 290 x 80 x 52 mm

Endoscope basket

**33-1062** 640 x 150 x 77 mm

Endoscope basket

**33-611** 274 x 172 x 43 mm

Wire basket

**32-993** 267 x 125 x 50 mm

Wire basket

**33-1070** 274 x 172 x 43 mm

Wire basket


**33-1073** 267 x 125 x 50 mm

Wire basket




Small fine mesh baskets with lid
Mesh distance = 2 mm x 2 mm

33-1100 30 x 30 x 30 mm**33-1101** 40 x 40 x 20 mm**33-1102** 40 x 40 x 40 mm**33-1103** 60 x 60 x 20 mm**33-1104** 60 x 60 x 40 mm**33-1105** 60 x 60 x 60 mm**33-1106** 80 x 80 x 20 mm**33-1107** 80 x 40 x 40 mm**33-1108** 80 x 80 x 40 mm**33-1109** 80 x 80 x 60 mm**33-1110** 80 x 80 x 80 mm**33-1111** 100 x 100 x 10 mm**33-1112** 120 x 120 x 30 mm**33-1113** 120 x 120 x 60 mm**33-1114** 120 x 120 x 80 mm**33-1115** 240 x 240 x 30 mm**33-1116** 240 x 240 x 60 mm**33-1117** 240 x 240 x 80 mm**33-1118** 240 x 240 x 100 mm




Specifications Sizes	h mm	 box and lid  box  lid		
1/1  540 x 253 x h	30	33-720	33-740	33-760
	50	33-721	33-741	
	70	33-722	33-742	
	100	33-723	33-743	
1/1  485 x 253 x h	30	33-724	33-744	33-761
	50	33-725	33-745	
	70	33-726	33-746	
	100	33-727	33-747	
3/4  405 x 253 x h	30	33-728	33-748	33-762
	50	33-729	33-749	
	70	33-730	33-750	
	100	33-731	33-751	
1/2  244 x 253 x h	30	33-732	33-752	33-763
	50	33-733	33-753	
	70	33-734	33-754	
	100	33-735	33-755	
Dental  275 x 175 x h	40	33-736	33-756	33-764
	95	33-737	33-757	
Small  265 x 125 x h	50	33-738	33-758	33-765



33-767
Stand
Height 20 mm



33-768
Stand
Height 14 mm



33-1119
Aluminium basket, with lid
245 x 130 x 55



33-1120
Aluminium basket, with lid
245 x 130 x 70



33-1121
Aluminium basket, with lid
245 x 130 x 105 mm



33-1080 600 x 400 x 150 mm

33-1081 600 x 400 x 300 mm

Laundry baskets



33-1085 600 x 400 x 150 mm

33-1086 600 x 400 x 300 mm

Laundry baskets


Drape	Sizes	for MODEL
	32-1600 130 x 60 cm 32-1601 140 x 100 cm 32-1602 185 x 45 cm 32-1603 180 x 60 cm	1/1
	32-1615 140 x 100 cm	3/4
	32-1630 90 x 60 cm 32-1631 90 x 80 cm 32-1632 130 x 45 cm 32-1633 130 x 60 cm	1/2
	32-1640 90 x 60 cm 32-1641 90 x 80 cm	Flat
	32-1650 60 x 40 cm	Mini
	32-1660 60 x 40 cm 32-1661 65 x 65 cm 32-1662 85 x 85 cm	Dental

OLUŞUM[®]



SURGICAL INSTRUMENTS CO., INC.

Sterile Container System



Accessories



The security plastic seals will be broken upon the lock is opened and this action makes it possible to follow up the sterilization process.



32-663
Security seal, made of plastic
100pcs. with indicator. (Blue)



32-662
Security seal, made of plastic
100pcs. with indicator. (Red)



32-664
Security seal, made of plastic
100pcs.(Blue)




32-660
Security seal, made of
plastic 100pcs.(Red)



32-661
Security seal, made of plastic
100pcs. with number.(Blue)




32-665
Security seal, made of plastic
100pcs. with number.(Red)



32-672
Labels, paper 100 pcs.
With indicator.
For full and half size containers




32-673
Labels, paper 100 pcs.
With indicator.
For small size and special
containers



- **32-700**
- **32-701**
- **32-702**
- **32-703**
- **32-704**
- **32-705**


Coloured identification labels
Engraved to your specifications
Maximum 16 letters



32-690
Dispenser for identification labels



32-680
Dispenser for identification labels




BLUE COLOURS

32-526 B	520 x 230 mm
32-528 B	470 x 230 mm
32-525 B	440 x 230 mm
32-524 B	380 x 230 mm

GREEN COLOURS

32-526 G	520 x 230 mm
32-528 G	470 x 230 mm
32-525 G	440 x 230 mm
32-524 G	380 x 230 mm




BLUE COLOURS

32-523 B	300 x 300 mm
32-527 B	230 x 220 mm
32-521 B	250 x 235 mm

GREEN COLOURS

32-523 G	300 x 300 mm
32-527 G	230 x 220 mm
32-521 G	250 x 235 mm



BLUE COLOURS

32-529 B	260 x 180 mm
32-520 B	275 x 125 mm


GREEN COLOURS

32-529 G	260 x 180 mm
32-520 G	275 x 125 mm




32-530	250 x 240 x 2.0 mm
32-531	480 x 250 x 2.0 mm

Silicone mat
Enables surgical instruments to be stored
without slip and dry well








34-100 L= 250 mm H= 41 mm
34-101 L= 152 mm H= 41 mm
34-102 L= 125 mm H= 41 mm
34-145 L= 97 mm H= 41 mm





34-103 L= 240 mm H= 41 mm
34-104 L= 150 mm H= 41 mm
34-105 L= 120 mm H= 41 mm
34-146 L= 90 mm H= 41 mm




34-106 L= 238 mm H= 41 mm
34-107 L= 150 mm H= 41 mm
34-108 L= 119 mm H= 41 mm
34-147 L= 87 mm H= 41 mm





34-109 L= 238 mm H= 41 mm
34-110 L= 158 mm H= 41 mm
34-111 L= 112 mm H= 41 mm
34-148 L= 77 mm H= 41 mm





34-112 L= 178 mm H= 20 mm
 Instruments holder, from Ø 3 to 8 mm




34-113 L= 73 mm H= 20 mm
 Instruments holder,
 Ø 5 mm





34-114 L= 73 mm H= 20 mm
 Instruments holder,
 Ø 5 mm




34-115 L= 178 mm H= 30 mm
 Instruments holder, from Ø 3 to 8 mm



34-116 L= 121 mm H= 18 mm




34-117 L= 123 mm H= 18 mm








34-124 L= 178 mm H= 20 mm
Instruments holder, from Ø 3 mm





34-125 L= 178 mm H= 20 mm
Instruments holder, from Ø 5 mm





34-123 L= 123 mm H= 18 mm





34-127 L= 23.5 mm H= 21 mm
Ø 8 mm






34-129 L= 23.5 mm H= 21 mm
Ø 5 mm




34-128 L= 23.5 mm H= 21 mm
Ø 13 mm




34-126 L= 23.5 mm H= 21 mm
Ø 10.5 mm







34-141 L= 72 mm H= 30 mm
Endoscope holder, double
one side Ø 15 mm open
one side from Ø 3 to Ø 15 mm




34-130 L= 72 mm H= 30 mm
Endoscope holder, double
one side Ø 15 mm open
one side from Ø 3 to Ø 15 mm




34-131 L= 72 mm H= 30 mm
Instruments holder
from Ø 3 to Ø 15 mm




34-132 L= 72 mm H= 30 mm
Instruments holder
from Ø 3 to Ø 8 mm




34-133 L= 38 mm H= 20 mm
Endoscope holder, double
one side 12 mm open
one side from Ø 3 to Ø 8 mm




34-134 L= 38 mm H= 30 mm
Endoscope holder, double
one side 12 mm open
one side from Ø 3 to Ø 8 mm




34-135 L= 38 mm H= 29 mm
Endoscope holder, Ø 4 mm closed




34-136 L= 38 mm H= 29 mm
Endoscope holder, Ø 6 mm closed




34-137 L= 38 mm H= 29 mm
Endoscope holder, Ø 8 mm closed



34-138 L= 38 mm H= 29 mm
Endoscope holder, Ø 12 mm closed



34-139 L= 38 mm H= 29 mm
Endoscope holder, Ø 12 mm closed



34-140 L= 38 mm H= 29 mm
Endoscope holder, Ø 15 mm closed



32-740

Single Textile
for half, middle
and full size containers




32-741

Single use paper filter 100 pcs.
With indicator. No marking



32-742

Teflon filter
Number of sterilization cycle > 1200



33-612

Single Textile



33-615

Single use paper filter 100 pcs.
With indicator. No marking



32-743

Teflon filter
Number of sterilization cycle > 1200



32-760

Single Textile



32-761

Single use paper filter 100 pcs.
With indicator. No marking



32-744

Teflon filter
Number of sterilization cycle > 1200




**32-720**

Filter retainer, for 1/1, 3/4, 1/2, XL
For disposable paper filter

32-721


For PTFE Teflon filter

**32-730**

Filter dispenser

**32-1715**

Filter Retainer, for dental model

**32-1708**

Filter dispenser

**32-1716**

Filter Retainer, for mini model

**32-1709**

Filter dispenser


SILICONE GASKETS FOR LIDS



	Model	Size
33-950	Mini	300 x 140 mm
33-951	Dental	310 x 190 mm
33-952	Flat	285 x 280 mm
33-953	1/2	285 x 280 mm
33-954	3/4	465 x 280 mm
33-955	1/1	580 x 280 mm
33-956	XL	590 x 370 mm


For mounting

- 32-710** Red
- 32-711** Yellow
- 32-712** Green
- 32-713** Blue
- 32-714** Grey
- 32-715** Black




Silicone handle has available in six different colors and detachable by user.



**33-1090**

Rack, for endoscopy
Holds 10 endoscopic instruments
Fits into a 1/1 size 20 cm deep container

Wire basket 1/1 50 mm
Container 1/1 200 mm
Silicone mat 520 x 230 mm

**33-1092**


Rack, for arthroscopy
Holds 20 arthroscopy punches
Fits into a 1/1 size 20 cm deep container

Wire basket 1/1 50 mm
Container 1/1 200 mm
Silicone mat 440 x 230 mm

**33-1094**

Rack, for Kerrison
Holds 10 Kerrison punches
20 cm deep container

Wire basket 1/2 50 mm
Wire basket 1/1 50 mm
Container 1/2 200 mm
Container 1/1 200 mm





32-716
Holder for identification labels
50 x 35 x 13 mm



32-717
Disposal stamp
40 x 35 x 9 mm
Requested number or letter
must be informed during
the order



32-691
Storage stand
390 x 85 x 153 mm




32-560 RED
32-561 BLUE
32-562 GREEN
32-563 YELLOW
32-564 GREY
32-565 BLACK



32-570
32-571
32-572
32-573
32-574
32-575

Coloured identification
labels Engraved to your
specification

Grey coloured recommended for use in Sterrad.
Color fading By effect of Hydrogen Peroxide Gas Plasma on other colors.



32-569
Information board with clip





32-600 40 X 5 X 15 mm
DIVIDER (use with **32-601** and **32-602**)



32-601 23 cm
32-602 47 cm
Sliding rails with securing screws



32-604 23 cm
32-605 47 cm
Limiting bar, with securing screws



32-607
Spiral holding device



32-609
Spiral holding device





32-615	11 mm
32-616	17 mm
32-617	23 mm
32-618	47 mm

Limiting bar




32-620
Spiral holding



32-622
Spiral holding



32-625
Spiral holding



Container System & Baskets

Classic & Bio-Barrier Model



OLUŞUM®

OLUŞUM CERRAHİ ALETLER SAN. VE TİC. LTD. ŞTİ.
OLUŞUM SURGICAL INSTRUMENTS INDUSTRY AND TRADE CO.
Aşağıçinik Beldesi Aşmek Sanayi Sitesi No: 5-6 Tekkeköy/SAMSUN
Tel : 90 (362) 264 81 97 - 264 83 17 Fax : 90 (362) 264 84 27
www.olusumcerrahi.com E-mail : bilgi@olusumcerrahi.com